




Infratonic S[®]

The Infratonic S[®] is excellent for emotional processing including:

- 1) Dissolving subliminal limiting beliefs*
- 2) Dissolving emotional cell memory*
- 3) Dissolving stuck emotions*
- 4) Dissolving stored emotional trauma in tissue*
- 5) Reaching chronic tension and inflammation stored deep in the body*



2

1

3

6

4

Infratonic S[®] Basic Operation

- 1. Charging** – When you need to charge the Infratonic S (I-S), and whenever the red LED ① blinks (indicating the battery is low), plug the charger into any socket (100-240V) and into the I-S at point ⑥ for 1.5 hours for a full charge, which will provide up to 50 hours of operation. The red LED will turn to yellow then green while charging and ~~will shut off automatically~~ when fully charged.
- 2. ON/OFF** – Hold down the middle button for 1 second to switch the I-S ON or OFF.
- 3. Timer** – To select the desired timer function, push the Left button repeatedly until you find the desired setting:
All Lights OFF = Continuous operation
Alarm Light ON = Beep will sound at the end of the 10 minute or 20 minute cycle.
- 4. Signals** – To select between the three signals, press the Right button until the desired light illuminates.

5. Basic Treatment Guidance – Apply the I-S to the area of congestion for 10 to 20 minutes (it is best to apply it for shorter periods of time when using initially). It can be applied directly to the skin, through clothing, or from a distance. The deep penetrating waves will also work through blankets, pillows, and plaster casts. The transducer can be used directly over metal pins and staples without risk of heating. If the signal feels too strong, just move the transducer a few inches to several feet away from your body until you feel comfortable. If you still feel discomfort of any sort, apply the Infratonic 9, Mobile Medic, or other Infratonic to soothe the emotional body.



Signal Description

Melt – This is most effective for dissolving subliminal limiting beliefs and other emotional sludge where there is no conscious connection to it.

Heart – Where there are heart issues associated with the residual emotional sludge to be dissolved, this signal is best.

Soothe – This signal is best for bringing buried emotional issues to conscious awareness, and is best where it is desired to learn the conscious patterns behind the emotions.

*Where you just want the discomfort to melt, or when you just don't want to get into where it comes from, start with **Melt**. Otherwise, start with one of the other two signals. It is often best to finish with **Melt** and to follow with the **Infratonic** or **Mobile Medic** for calming and smoothing.*

x

Treatment Precautions

- Do not use directly against extremely inflamed or infected skin.
- Do not use directly on the leg calf if blood clots are suspected.
- Keep the I-S at least 6 inches away from implanted electronic devices, including pacemakers.
- Do not use the I-S while operating a motor vehicle or heavy equipment.

Specifications

Treatment	Infratonic massager
Frequency	8-800Hz
Battery Specs	3.6V, 700mAh NiCd rechargeable battery
Battery Life	Up to 50 hours
Charger	100-240V, 50/60Hz, 9V
Charging Time	1.5 hours for full charge, automatically shuts off when charged
Recharge	When BAT light is lit



*27130A Paseo Espada, Suite 1407
San Juan Capistrano, CA 92675*

(949) 361-3976

(800) 743-5608

www.soundvitality.com