



The Mobile Medic®

Portable. Effective. Convenient.

With the lightweight, battery powered **Mobile Medic®**, you can enjoy all the benefits of Infrasonic therapy anywhere. Be sure to carry it if you have any condition with recurring pain or discomfort. Also take it along to sporting events and during travel, to provide early treatment of unexpected injuries and illnesses. Its deep, penetrating, soothing waves can be applied directly to the skin or through blankets, pillows and plaster casts, and can be gently moved between locations.

Using the Mobile Medic®

- 1. Charging:** When you first receive the ***Mobile Medic®***, and whenever it needs charging, plug the charger into any wall socket (100-240V, 50-60Hz) and into the ***Mobile Medic®*** for 16 hours for a full charge. You will get about 40 hours of operation per charge.
- 2. Operation:** Simply switch the switch for ON and you will feel it begin to vibrate. If it doesn't vibrate or if it makes a high pitch sound, it needs to be charged. When treatment is done, switch the switch OFF.
- 3. Application:** Apply the ***Mobile Medic®*** to the pain or discomfort for 10-20 minutes, or for as long as you like to provide extraordinary pain relief in cases where nothing else has worked before, and to provide emergency relief when nothing else is available. Apply to the thymus (over upper sternum) for deep relaxation and tranquil sleep. Don't worry if it runs all night. It can run every night for a week before it needs to be recharged!

4. Treatment: Simply put the **Mobile Medic®** where it hurts, when it hurts – including directly on the skin (unless it is extremely inflamed or infected). The deep penetrating waves will also work through blankets, pillows and plaster casts. The transducer can be used directly over metal pins and staples without risk of heating. If the signal feels too strong, just move the transducer a few inches to several feet away from your body until you feel a gentle stimulating effect.

Treatment Precautions

- *Do not use directly against extremely inflamed or infected skin.*
- *Do not use directly on the leg calf if blood clots are suspected.*
- *Keep the transducer at least six inches away from an implanted device, such as a pacemaker or defibrillator.*

Specifications:

Treatment: Infratonic Massager

Frequency: 8.5 – 13.8 Hz.

Battery Life: Up to 40 hours

Charger: 100-240V 50-60Hz 9V



27130A Paseo Espada, Suite 1407

San Juan Capistrano, CA 92675

(949) 361-3976

www.soundvitality.com

Copyright 2011. Sound Vitality. All Rights Reserved.