

INFRATONIC 8000	VITALIZE SUBCELLULAR HEALING	RELEASE CELLULAR TRAUMA	EVOKE	RELAX TIGHT MUSCLES	REDUCE STAGNANT ALPHA	UPLIFT	REDUCE PHYSICAL PAIN	REDUCE WORRY
SLEEP	YES	YES		YES	YES	YES	YES	
Pain		YES	YES	YES	YES	YES	YES	YES
FULL	YES	YES	YES	YES	YES	YES	YES	YES
INFRATONIC 8				New NanoSynthesis Signals				
Body	YES	YES						
MIND		YES						
SPIRIT	YES	YES	YES					

The Infratonic 8000 combines eight separate signals

Infratonic 8000

New for 2006 with 8 signals combined through NanoSynthesis technology, the Infratonic 8000 accelerates recovery like never before!

Infratonic Technology.... Pg 2

8000 Upgrade Offer......Pg 3

Product Line-up.....Pg 3

Immune Enhancement..Pg 4

The NEW Infratonic 8000 was designed to maximize accelerated recovery and Mind/Body healing. This was made possible because we determined that we have two distinct markets, Pain Management, and accelerated recovery + mind/body healing.

The SoniCalm 3 (released in January) optimizes the pain management effects of the Infratonic technol-

INSIDE:

ogy. Now, we have released Infratonic 8000 for accelerated recovery and increased vitality. While the 8000 is ef-

fective at acute and chronic pain, it does have some accelerated recovery

signals that can increase sensations of restlessness and pain, thus it is not strictly a pain management device.

The new technology that makes the 8000 possible is presented on page 2. It involves new ways to generate and combine unpredictable signals for maximum penetration and effectiveness. Eight distinct NanoSynthesis signals are combined in the 8000

> through the process of Sonic Sequencing, which makes each signal more powerful than it would be alone by moving

traumatized tissue, gently but unexpectedly, from one vibrational state

to another. A new chaotic carrier is applied to these signals that makes them more unpredictable and, in addition, is a therapeutic signal to remove mind/body barriers to healing and to lift mood.

The new 8000 also simplifies application of Infratonic Therapy by providing three powerful settings, SLEEP, PAIN and FULL. For best results in accelerated recovery, use the FULL setting. It combines all of the NanoSynthesis technology, not just the three signals in the Infratonic 8 (Sub-cellular Vitalization, Calm Cellular Trauma, and Insight) but also signals to relax tight muscles, uplift intention, reduce physical pain, and reduce worry.

FULL is a powerful setting and may result in some discomfort. If you start feeling restless or antsy, switch to the SLEEP setting. If you start feeling increased pain, switch to the PAIN setting and use the FULL setting a little more each day. You can also hold the transducer several inches or feet away from the body, or sweep with the transducer out the hands and feet.

Upgrade now from Infratonic 8 to the 8000. See page 3 for details.



What is so Unique about CHI's NEW Infratonic Therapy?



Sound affects physical structures


Sound Therapy is the use of sound vibrations for a positive affect on human health. Familiar to most are the therapeutic affects of music – recognized across many cultures for its ability to reduce stress and inspire emotion.

Ultrasound Waves are high frequency sound waves above the frequency of 20,000 hertz. Imperceptible to the human ears, ultrasound therapy provides diagnostic imagery in healthcare, and is often used in joint and muscle therapy.



Whales communicate with infrasound


Infrasound Waves are ultra low frequency sound waves with frequencies below human perception below 30 hertz. Unlike ultrasound waves, infrasound waves travel great distances without dissipation by physical obstacles. These are the frequencies that whales use to communicate, that scientists study to predict earthquakes, and that neurologists read on EEG in the study of brain waves and human thought.



The Chaos Therapy Signal


Infratonic Therapy[™] uses specific, randomly generated bands of infrasound waves to penetrate superficial and deep layers of the body. It works by interrupting trauma at the cellular level, using the same frequencies that the body itself uses for intercommunication. Infratonic therapy's ultra-low frequency sound waves increase micro-circulation, reduce swelling, relieve chronic pain, and promote healing.

Chaos Therapy® -- Introduced in 1997, this patented method infuses the body with highly unpredictable signals of specific bandwidth for greater penetration and higher effectiveness.



The CHI NanoSynthesis computer produces Chirp Entrainment signals

- 1. NanoSynthesis™— The first of four patent pending Infratonic technologies introduced in 2006 generates a highly variable therapeutic signal point-by-point with continuous phase angle and abrupt jumps from frequency to frequency within the desired band. The others are:
- 2. Chirp Entrainment[™]-- Short chirps of rising frequency are inserted into a higher frequency signal to entrain lower frequency signals and lift them into the desired range. Used to lift pain signals from the physical to the emotional range, and to lift worry signals from the analytical into the intuitive range.



The FULL setting of the 8000 combines eight Infratonic Signals

- 3. Sonic Sequencing[™] -- Switching seamlessly from signal to signal, repeatedly repositioning the body's vibrational system for maximum effectiveness.
- 4. Chaotic Carrier[™] -- Unpredictable high frequency signal that rides on the NanoSynthesis signal to increase its unpredictability. This is also the first Infratonic signal that addresses intention, lifting the frequency of the body's vibrational field to reduce limitations to accelerated healing.

PRODUCTS: | INFRATONIC 8000:

Palm Transducer 3:

The long anticipated 3rd Generation Palm Transducer is almost here. We expect to release it sometime in April 2006. With Linear Pulse Modulation and a strong solid feel on the body it is better than ever. It fits comfortably in the palm of your hand, has a slim profile that allows it to fit comfortably under you when

on a bed or a comfortable chair, and its grabby rubber feet and low center of gravity hold it firmly in place on a



patient's back. Black or White. Supplies are limited.

Palm Transducer: \$395

CHIonizer:

The CHIonizer is an Air Travel Vitalizer that keeps you vitalized and comfortable on long flights and kills airborne pathogens with a gentle flow of ionized air and ozone. Enjoy your flight and stay healthy!! Use it wherever you are faced with low vitality or people who might be coughing and sneezing: offices, buses, and polluted environments.



Great for shopping malls as well. Don't leave Earth without it!! CHIonizer: \$119

The 8000 is is our most advanced accelerated recovery



device, It combines 8 unique Infratonic signals for accelerated recovery on a wide variety of levels including signals that provide increased sub-cellular vitality, reduced physical pain, cellular cooperation, reduced cellular trauma, reduced worry, increased insight, and lifting of intention. These signals combine to make this the most effective therapeutic massager on the planet.

Infratonic 8000: \$895

SoniCalm 3:

Extraordinary pain relief is what people across the country are reporting with the SoniCalm 3. To accommodate the separate needs of Pain Management and Accelerated Recovery, CHI developed two separate units. The SoniCalm 3 excludes those Infratonic signals that vitalize subcellular activity and dissolve



mind/body barriers to healing, but that can increase discomfort. So where pain is the primary indication, or where the doctor's practice specializes in pain, the SoniCalm 3 is the best answer.

SoniCalm 3: \$695

Upgrade to the Infratonic 8000:

To upgrade to the Infratonic 8000, just send in your Infratonic 8 or Infrasound 8 with \$189.

Special, Limited Time Special Offer #1:

If you purchase our new palm transducer as a part of your upgrade, you will save \$200 from the regular Palm transducer price of \$395.

Special, Limited Time Special Offer #2:

If you purchase CHIonizer as a part of your upgrade, you will save \$20 off the regular price of \$119.

To upgrade to the Infratonic 8000, just call us today at

1 (800) 743-5608

Offer expires May 31, 2006

A Strong Immune System: Quick Fixes and Long-Term Answers

With long lasting strains of flu floating around that can cause missed weeks from work, and threats of newer and more potent epidemics like the Avian Flu lurking ominously around

the corner, we need to be prepared. We can stock up on flu medicines, but that's a weak answer. Diet, exercise, and nutritional supplements all boost immune competency in important ways, but many of us need more.


What we need is a strong, competent immune system that will quickly identify the virus and build an army of antibodies to defeat it quickly. Three

approaches are described here. First, we can restore immune competence by reducing trauma in the immune system. Second, we can attack airborne viruses before they get to us, and third, we can augment the electrical vitality the immune system needs to function strongly in a protracted battle against the flu.

Immune Enhancement Protocol:

Because fear and anxiety cause immune competency to plummet, the objective of this protocol is to detraumatize the immune system so it will be at maximum effectiveness when it is needed. Simply apply the Infratonic for ten minutes three times a week to the center of the chest just

below the collarbone, and relax. The FULL signal of the Infratonic 8000 is the best, though all Infratonic signals will help. If you use it just before bed, use the SLEEP signal in-



Blood serum antibody level of rabbits shows accelerated immune response in the group treated with the Infratonic.

approaches are described here. First, stead. Our goal is a reduction from 10 we can restore immune competence by reducing trauma in the immune system. Second, we can attack air
stead. Our goal is a reduction from 10 to 3 sick days from flu per year. Please let us know if this protocol is helpful for you!!

Copyright © 2006 CHI Institute. All Rights Reserved.

For an immune system boost that can make the difference between catching and fighting off a virus, the CHIonizer can be worn whenever you are in an enclosed environment where airborne pathogens (sneezy fellow passengers, shoppers, or coworkers) are present.

The CHIonizer infuses the body with electrical vitality both by producing ionized air that vitalitzes the body and kills airborne pathogens directly, and by providing an ionic discharge within the body's electrical field. (within about ½ inch from the skin)

When you have the flu, the CHIonizer is also very valuable to wear. I had a chronic cough for three weeks, then got on a plane for Australia. I put on the CHIonizer as the plane took off, and coughed for the first two hours. Then it stopped. I had no coughing for the rest of my trip!!! To test this out, I found five people who had suffered from chronic cough for several days. For every one, the cough disappeared within two days of starting to wear the device. The CHIonizer is clearly valuable where the immune system is run-down or overloaded.

See page 3 for CHIonizer Details



Presorted Std. U.S. Postage PAID Santa Ana, CA Permit No. 1135